

TQSoft and TQAero

Software for the 1586A Super-DAQ or 2638A Hydra Series III used in thermal validation studies to meet the data collection and reporting requirements of regulated industries

TQSoft – For pharmaceutical / biotechnology companies:

- United States Food and Drug Administration (FDA) **21 CFR Part 11** regulations on electronic records and signatures for incubation, sterilization, freezing, drying and temperature mapping validation applications
- European standards for sterilization, decontamination, and disinfecting (EN554, EN285, EN15883, HTM2010, HTM2030)
- ISO 15833 requirements for washer-disinfectors and ISO 17025 competence requirements for testing and calibration laboratories

TQAero – For heat treatment service / aerospace / auto / transportation companies. Complies with SAE International **AMS 2750** requirements for metal heat treating applications



1586A Super-DAQ with DAQ-STAQ Multiplexer



2638A Hydra Series III



21 CFR Part 11



TQSoft: Thermal study example

1. Set-up equipment


installed on PC or laptop



Temperature source

- 9190A Field Metrology Well

Reference probe

- 5615 Secondary Reference PRT

Test sensors

- Type T thermocouples



Data acquisition system

- 1586A Super-DAQ with DAQ-STAQ

Other Fluke equipment supported:

Data acquisition system

- 2638A Hydra Series III and previous generations

Thermometer readouts

- 1502A/1504 Tweener
- 1523/1524 Reference Thermometers
- 1560 Black Stack

Temperature sources

- 9142/9143/9144 Field Metrology Wells
- 9170/9171/9172/9173 Metrology Wells
- 6102/7102/7103 Micro-Baths
- 6330/7320/7340/7380 Compact Baths
- 6331/7321/7341/7381 Deep-Well Compact Baths
- 9150 Thermocouple Furnace
- Other drywell calibrators

Additional equipment supported from Agilent, Anville, Eurotherm, Kaye, TMI, Yokogawa, Ametek, and Isotech

2. Automatically calibrate sensors

- Enter stability criteria
- TQSoft communicates with temperature source, reference probe, and test sensors
- When reference and sensors are stable as defined, TQSoft automatically applies all correction offsets
- It also generates and archives a full and traceable calibration report




Calibration Progress dialog box showing a table of stability requirements and a progress bar. The table has columns: Channel No, Value (°C), Countdown on Stability Requirements, Biggest Drift last 60 seconds, and Deviation from Reference. The table data is as follows:

Channel No	Value (°C)	Countdown on Stability Requirements	Biggest Drift last 60 seconds	Deviation from Reference
1	45.1	03:00	0.4	0.7
3	45.8			1.4
-	44.33	00:46	0.40	

3. Set-up thermal study for chamber



Stability chamber for pharma/biotech products



- Use TQSoft to manage a 3D wire diagram of the chamber space to be studied
- Position sensors in the correct locations
- Rotate and magnify position diagram. Edit if needed.
- Sample test specifications are included that can be used while learning TQSoft

4. Test and view results



installed on PC or laptop



Stability chamber with Test sensors





Data acquisition system
• 1586A Super-DAQ with DAQ-STAQ

- Type T thermocouples

- Temperature measurements can be displayed during testing
- Temperature gradients with color shading at sensor positions can be shown
- After testing, replay with a time slider showing temperature behavior in the chamber

TQSoft: Thermal study example cont.

5. Create report



- TQSoft includes a “Report Manager” for design and printing of reports
- Includes a large library of report templates for all processes
- Reports can include graphics, complex calculations based on acquired test data, lethality formulae, and engineering lookup tables. The calculations and lookup tables can be validated.
- Reports are presentation quality, quickly produced and secure

Instrumentation Check During Exposure Time				
Drain Temperature				
Indicated	Recorded	Measured	Difference	Limit
Hold Indic.	Hold Reco	Hold Meas	Hold MRF	Hold MRF
Chamber Pressure				
Indicated	Recorded	Measured	Difference	Limit
Hold Indic.	Hold Reco	Hold Meas	Hold MRF	Hold MRF
Timer Accuracy				
Exposure Time				
Indicated	Recorded	Measured	Difference	Limit
Hold Indic.	Hold Reco	Hold Meas	Hold MRF	Hold MRF

Cycle Start Time	Cycle Start Time
End of Cycle	Cycle end time
Cycle Time	Total cycle duration
Start Hold Time	Sterilising End of Equ
End Hold Time	Sterilising end time
Hold Time	Sterilising hold time
Start Exposure Time	Sterilising Start of Equ
End Exposure Time	Sterilising end time
Exposure Time	Sterilising plateau time

TQSoft and TQAero Ordering Information



Model Number	Description	What's included
TQSOFT-IQ/OQ	TQSoft Pharma Process Analysis and Validation Software	Data acquisition software with activation key code. Electronic signatures and operator action audit trail. Configurable for use in discrete process validation. Provides FDA 21 CFR part 11 compliant user access controls. Includes IQ/OQ Documentation to help customers test their validation system with their processes. Email support and one remote installation and training session is included.
TQAERO	TQAero Furnace Validation Software	Data acquisition software with activation key code. NADCAP compliant user access controls and operator action audit trail. Configurable for use in discrete process validation. TQAero is a configuration of TQSoft providing complete support for AMS2750 requirements. Email support and one remote installation and training session is included.
VAL-BNDR-TQS	Validation Reference Binder	A set of TQSolutions documents on CD summarizing TQSoft system acceptance testing to GAMP standards.